

Synchronization Improves Potentiostat Performance

Robert S. Rodgers and Gregory A. Martinchek

Gamry Instruments, Inc.

Warminster, PA



Noise in Potentiostat Systems

- Intrinsic Noise
 - Generated by electrochemical cell
 - Electrochemical Noise
- External Pickup
 - Mains frequency and harmonics
 - 50/60 Hz, 100/120 Hz
 - Nearby RF
 - Motors, CRTs, other instrumentation
- Move the lab?
- Shielding
 - Faraday Cage
 - Surround cell and leads with conducting and *earthed* shield, connect to ground of instrument

Noise in Potentiostat Systems

- Internal Noise
 - Generated within the Potentiostat system
 - Component noise
 - Operational Amplifiers (Op-Amps)
 - Resistors (Johnson noise)
 - Voltage reference (Internal standard)
 - Power supply
 - Often a major contributor
 - Limits smallest signal that is measurable
 - pA – aA currents, μV voltages
 - Applied Potential Noise
 - Actually applied to the cell
 - Noise can hide natural effects – EIS, ECN
 - » Very low amplitude signals
 - May yield large currents on capacitors (double layer!)

Internally Generated Noise


- Power Supply
 - Linear
 - “Good old days”
 - Most Noise & ripple @ mains, 2x, 3x mains
 - Internally radiated
 - Current spikes on ground plane
 - Controlling the noise
 - Filtering
 - Also filters **signal** !
 - Integrating over a mains cycle
 - Also filters **signal** !
 - Take data slowly!

Internally Generated Noise

- Power Supply
 - Switch-Mode or Switching
 - High efficiency (less heat)
 - Can serve as DC-to-DC converter
 - Noise occurs at switching frequency and harmonics
 - YOU can control the frequency!
 - Above frequencies of interest
 - Controlling the noise
 - **Synchronize** data acquisition/waveform generation to the switching frequency

Switch-Mode Power Supply


- Input voltage can be DC
 - DC-to-DC converter
- May provide isolation between Input and Output:
 - Special Transformer
 - Control circuitry for Switch S must also be isolated
 - Opto-coupler
- Isolation allows operation with **grounded/earthed** electrodes!



Switch-Mode Power Supply as used by Gamry

Overall Schematic


- Isolation of Potentiostat provided by
 - Digital USB isolators (USB may be earthed)
 - Power Supply Transformer
 - Opto-isolator on PS Controller Clock
 - Experiments with **Earthed electrodes** are possible



Isolation in the Reference™ 600

One Crystal Controls All

- A single 32 MHz crystal is the master timing for nearly everything!
- Power supply switching frequency
 - Opto-isolated to PS controller
- Data acquisition *via* "State Machine"
 - Synchronize with internally generated noise
 - Always sample at the same instant of internal noise waveform during an experiment



Clock generation in a Reference™ 600


State Machine

- **Finite state automaton**
 - “state”
 - Has outputs
 - Digital signals to control data acquisition
 - Has inputs
 - Digital signals from the external world, or even its own outputs!
 - Has rules for transitions into another state
 - Examples
 - Door: Open or Closed
 - Traffic Signal: Green, Yellow, Red, Yellow-Red




State Machine

- Implemented as an integrated circuit
- A CLOCK synchronizes all transitions
 - Prevents ambiguities during transitions
 - Controls time between transitions/decisions
- Controls data acquisition timing
 - Within a point
 - Current interrupt
 - When to take data
 - Time between points




State Machine

- Because power supply noise frequency (300 kHz) and data acquisition are derived from same clock, they can be made synchronous
- Data acquisition
 - 3,333 μs or a multiple!
- Always measure data at the same point on (possible) noise waveform



Noise Measurements

- Reference™ 600 has two “synchronization” digital signals
 - Input pulse to start data acquisition
 - Pulse out when data acquisition starts
- Effect of Non-synchronized?
 - Two experiments, same data point timing
 - (1) Normal clocking - use **A** only
 - (2) **A** is triggered by pulse from **another** Reference 600 (**B**)



Noise Measurements


- Measured V across 2000 ohm resistor @ 0V with potentiostat **A**
 - Faraday Cage!
 - Same potentiostat, but trigger came from A or B
- 2000 data points averaged

Point Timing/ μ s	Internal Trigger (A) (synchronized)		External Trigger (B) (un-synchronized)	
	Average/ μ V	Std Dev/ μ V	Average/ μ V	Std Dev/ μ V
10,00	35,1	69,6	-43,6	124,1
"	50,9	68,8	0,9	123,7
3,33	426	80,9	-131	370
"	348	80,3	-168	383
3,33 *	262	47,6	528	92,2
" *	261	47,1	449	96.8

* 200 kHz filter applied


Synchronization

- Synchronization of power supply and data acquisition improves measured noise by a factor of 2 or more
- 200 kHz filter did not entirely remove noise with external trigger
 - Noise >200kHz is attenuated, reduces noise bandwidth
 - Synchronization does not help component noise
 - Possibly noise was introduced post-filter




Synchronization and EIS

- Conventional Wisdom
 - You **MUST** sample many times/cycle
 - You **MUST** use an “anti-Aliasing” filter to eliminate all frequencies above the “Nyquist” Frequency
 - Nyquist Frequency = $f_{\text{Acquisition}}/2$
- If you wish to measure EIS @ 1 MHz, you must take data every 500 ns or **faster!**
- Expensive: € € € € €



EIS by Sub-Harmonic Sampling


- Conventional Wisdom
 - Assumes you know **nothing** about the signal
- EIS
 - Apply a signal of **known** frequency
 - Measure E, I at **that frequency**
- Break the “rules”
- Can use aliasing to your advantage!
 - **Only** if data acquisition is synchronized with waveform generation!



One **data point** every $\sim 1,05$ cycles of **sine wave**


EIS by Sub-Harmonic Sampling

- Traditional waveform circuits are unsuitable
 - Accuracy controlled by R, C, L values
 - Inaccurate
 - Frequency can drift
- DDS – Direct Digital Synthesis
 - Produces a low distortion analog sine wave
 - Programmable frequency
 - $24\text{MHz} * N / 2^{32}$
 - Used in TV, radio, cell phones




What is Sub-Harmonic Sampling?

- Requires a repetitive signal
 - Box-Car Averaging
 - Down-mixing
 - Heterodyne
- Example: $f = 2\,777,777\dots$ Hz; period = $1/f = 360\ \mu\text{s}$
- Take a data point every $396\ \mu\text{s}$ ($f / 1,1$)
 - < 1 point per cycle
 - $f_{\text{ACQUISITION}} >$ Nyquist frequency
 - Point 1 396° (36°); point 2 792° (72°); point 3 1188° (108°); ... point 10 3960° (360° or 0°)
- After $3960\ \mu\text{s}$ we have 10 points, evenly spaced every 36°



Sub-Harmonic Sampling

- At higher frequencies, **many** cycles may be skipped
 - Start skipping cycles above ~ 10 kHz
 - @ 1 MHz, skip 166 cycles (of 1 MHz) between readings
- Only possible if DDS and data acquisition clock are locked together
 - @ 1MHz: $10\text{ns} \rightarrow 3,6^\circ$
- Stringent requirements on Sampling ADC as well
 - **ADC input** should have 1 MHz bandwidth!
 - Not often specified!
 - Sample-and-Hold “Aperture” jitter is an important specification
 - 5 ps rms typical (Reference™ 600)





Is Internal Clock Necessary?

- Two Reference™ 600 Potentiostats
 - Normal connection
 - Internal DDS of **A**
 - External DDS
 - Same DDS type
 - Same crystal type
 - Signal Out from **B** → External Signal In **A**
 - Measured EIS 65 Hz-357 Hz
 - Applied Vac quite similar for both experiments



Is Internal Clock Necessary?

	$\diamond(Z)$	$\diamond(\angle)$
Internal	0,007%	0,004°
External	0,14 %	0,072°



YES!

Acknowledgments

- Greg Martinchek,
 - Co-author and hardware architect @ Gamry
- Max Yaffe
 - Software architect @ Gamry and architect of Sub-harmonic sampling
- The hardware and software engineers at Gamry
 - For tolerating my thousands of questions!
- You, the audience
 - Danke sehr für Ihre Aufmerksamkeit !!